




**3 POSITION GERICHAIR  
TNR630  
OWNERS MANUAL**



**Seat Width: 21”  
Weight Capacity: 250 lbs.**



**ATTENTION**

Before using this Medical Device read this manual. If you are unable to understand, contact your equipment provider for technical support before attempting to use this product.

## **SAFETY PRECAUTIONS**


### **CAUTION:**

Do not operate recliner by applying force or pressure to the leg rest. Avoid placing hands or legs in creases near recliner mechanism.

- ⚠ Periodically inspect all parts and replace as needed.
- ⚠ Occupant weight should not exceed 250 lbs.
- ⚠ Caster wheels should be locked during occupant transfer.

## **PRODUCT DIAGRAM**

### **3 POSITION GERICHAIR TNR630**



- |                                     |   |
|-------------------------------------|---|
| <b>A.</b> Upholstered Back          | <b>B.</b> Back Assembly                         |
| <b>C.</b> PVC Push Bar              | <b>D.</b> Holes in Backrest                     |
| <b>E.</b> Holes with Screws         | <b>F.</b> Recliner Lock Bar Storage             |
| <b>G.</b> Table Storage             | <b>H.</b> Lock-Release Pull Button              |
| <b>I.</b> Holes in Horizontal Brace | <b>J.</b> Label Specifying Date of Manufacturer |
| <b>K.</b> Locking Pins              | <b>L.</b> Lower Rear Support Tube               |

## **ASSEMBLY**

1. Remove back assembly from unassembled chair and aside momentarily.
2. Slide the 1 piece back assembly into the square tubes of the lower frame assembly so that the holes in the back rest line up with the holes located on horizontal brace of the lower frame assembly.
3. Raise recline lock bar adjust position of back assembly until lock release button positively engages middle slot lock bar link when chair is in storage position.
4. Insert screws through holes of the horizontal brace and tighten all screws with screwdriver.

**NOTE:** Make sure that the top screw on both sides of the coated square tubing from the back assembly are sufficiently secure before using recliner.

## **OPERATION**

### **PLACING THE UNIT IN RECLINING POSITION BY THE OCCUPANT**

1. Place hands on armrests and push against seat back. Leg rest will automatically adjust to each recline position.

**NOTE:** The table can only be used in upright sitting position. Patient should be transported in reclined position with recliner lock set.

### **PLACING THE UNIT IN RECLINING POSITION BY THE ATTENDANT**

1. Ensure that recliner lock bar is in stored position.
2. Stand at either side of the recliner.
3. Grasp armrest with one hand and PVC push bar.
4. Push down on push bar to recline the back to the first recline position. To achieve second recline position push on push bar again.
5. To raise upright, pull forward on push bar.
6. Set recliner Lock Bar.

## RECLINING LOCK OPERATION

### To Lock:

1. Place recliner into desired position.
2. Pull out Lock-Release Pull Button to release Lock Bar from stored position.
3. Rotate recliner Lock Bar down until notches are fully engaged over Locking Pins.

### To Unlock:

1. Lift recliner Lock Bar up to disengage from Locking Pins.
2. While holding Lock-Release Pull Button out, rotate recline Lock Bar all the way up until completely against seat back.
3. Release Lock-Release Pull Button to hold recliner Lock Bar in stored position.


**NOTE:** Do not attempt to change position of chair with recliner lock bar engaged. This could cause damage to mechanism.

## OPERATING THE BLOW MOLDED TABLE

### WARNING:

Casters should be in locked position during patient transfer.

1. Underneath each arm, pull down and hold pull button and slide Lockout Keys forward on both sides.
2. Pick up table by circular table arm and slide into receiving tubes in front of each arm.
3. Table can slide on and off by keeping lock/lockout key forward.
4. Table can be locked into one of five positions. Pull down the pull button, hold the pull button and slide the key backwards. Proceed to the other arm and move that key backwards. Table will now be locked.



Articulating headrest may be adjusted to desired position by pulling forward or pushing back.

**STORAGE:** The blow molded table stores in small round tube underneath the left armrest, leaving the tabletop facing out.

# **MAINTENANCE**

## **CLEANING THE CHAIR (UPHOLSTERY MAINTENANCE)**

- Ordinary dirt: Wash with warm water and a mild soap or detergent.
- Stubborn Dirt: Using a soft bristle brush, follow directions above.
- Ground-in Dirt: Scrub with powdered cleanser or similar detergent.
- Chewing gum: Scrape carefully and remove with kerosene or naphtha.
- Nail Polish and Remover: Blot immediately to prevent damage.
- Tars, Asphalt, Creosote: Remove to prevent staining; Clean area with kerosene or naphtha.
- Paint, Shoe Polish, Heel Marks: Remove immediately; Use white cloth dampened in kerosene, naphtha or turpentine. Do not use paint remover or liquid brush cleaners.
- Ball Point Pen: Blot immediately with a white cloth dampened in water or mineral oil.
- Washing and Refinishing: Improve wearability. Use any quality paste wax. A diluted solution of 10% (or less) of bleach or Clorox may be used when cleaning this chair.

## **DETECTING WEAR AND DAMAGE**

- Check every 3 months that all nuts and screws are tight and have not loosened.
- Check and clean every 3 months the casters and the brakes on the casters.
- Check and clean every 3 months the Lock Recline Mechanism.
- Contact your equipment provider for replacement of any defective parts immediately and ensure that the chair is not used until repair job is completed.

### **Keep For Your Records**

Name of dealer where you purchased:

---

Dealer's phone number where you purchased

---

Serial number ( if applicable ):

---

Please keep this information for your records to help us in assisting you if you need a part replaced.

**Thank you.**

**West Coast Distribution**

**Tel 800-367-6160**

**East Coast Distribution**

**Tel 800-548-6596**

## LIMITED WARRANTY

The warranty below has been drafted to comply with Federal Law applicable to products manufactured after July 4, 1975. This warranty is extended only to the original purchaser/consumer (or dealer non-consumer who does not buy for resale).

Tuffcare® warrants that its proprietary chair will be free from defective workmanship and materials for a period of one (1) years following the date of original purchasing date with the following exception: Parts on the chair are warranted for six (6) months. The warranty period commences on the original purchasing date. If within such warranty period any such product proven to Tuffcare® satisfaction to be defective, such product will be repaired or replaced at Tuffcare® option.

Tuffcare® sole obligation and your exclusive remedy under this warranty shall be limited to such repair and/or replacement. This warranty does not include any labor charges incurred in replacement parts installation. Freight charges to the factory are at the expense of the consumer or seller. Return freight charges will be prepaid by Tuffcare®. For warranty service, please contact the authorized dealer from whom you purchased your Tuffcare® products. In the event that you do not receive satisfactory warranty service, please write directly to: Tuffcare®: 4977 La Palma Ave, Anaheim, California, 92807. DO NOT RETURN PRODUCTS WITHOUT PRIOR AUTHORIZATION.

**LIMITATIONS AND EXCLUSIONS:** The foregoing warranty shall not apply to products subjected to negligence, abuse, misuse, improper operation, improper maintenance, improper storage or damages beyond Tuffcare® control. The evaluation will be solely determined by Tuffcare®. The warranty shall not apply to problems arising from normal wear, or failure to follow instructions, or if parts are not manufactured by Tuffcare®, or do not comply with original equipment specifications that are added to Tuffcare® product.

THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE AND SHALL NOT EXTEND BEYOND THE DURATION OF THE EXPRESS WARRANTY PROVIDED HEREIN. TUFFCARE® SHALL NOT BE LIABLE FOR ANY CONSEQUENTIAL OR INCIDENTAL DAMAGES WHATSOEVER.

This warranty gives you specific rights and you may also have other legal rights which vary from state to state. Some states do not allow the exclusion or limitations of incidental or consequential damage, or limitation on how long an implied warranty

**NOTES:**

[illegible]