Comfy - Aire AAP501

UPHOLSTERY COLOR: NAVY BLUE

Mattress Base: BLACK



Our AAP501 alternating overlay system is designed to fit over an existing mattress to deliver pressure reduction. This product offers "at risk" patients the most comfortable support surface is its class and at the same time provides an environment that reduces the chance of developing pressure ulcers. Key features include 24 alternating air cells, as well as a CPR quick release valve for comfort and safety.



Comfy - Aire AAP501

UPHOLSTERY COLOR: NAVY BLUE

Mattress Base: BLACK

STANDARD FEATURES

- * Fits over existing mattress
- * CPC connectors
- * CPR quick release Valves
- * Zipper top cover for easy removal
- * Vapor permeable top cover which is machine washable
- * Vinyl covering on air tube
- * LED indicator
- * Weight limit: 350 lbs


OPTIONAL ACCESSORIES

* Carrying bag

* Fast filling air pump

SPECIFICATIONS

PUMP SPECIFICATIONS:		MATTRESS SPECIFICATION:	
Power:	110V/ 60Hz., 230V/ 50Hz., 5.2W	Cell Depth:	5 1/4" Inflated
Pressure Range:	30 - 50 mmHg.	Air Cell:	20 Alternating Air Cells
Air Volume Output:	5 liters on 110V/ 6 liters on 230V	Top Cover:	Removable & Washable
Cycle Time:	5 mins. (50Hz & 60Hz)	Foam Base:	None
Power Cord:	Hospital Grade Cord & Plug	Odor Control Cover:	Yes
Power Cord Length	12 feet	Air Tube Covering	Yes
Firmness Adjustment:	Yes	CPR Deflation:	Yes
CPC Connectors:	Yes	CPC Connectors:	Yes
Dimension:	8 1/2" x 8 1/2" x 4"	Inflated Dimension:) " x 35" x 80"
Weight:	5.3 lbs.	Weight:	14 lbs.
Agency Approval:	UL, CSA		



Inflation and deflation of air cells in alternating cycles provide pressure relief for patient's sensitive tissues.



Total of 2€ 5 F/4" air cells whichÁ provide the necessary lift toÁ achieve pressure relief.



High capacity 5 liter air pump is the quietest on the market. Features an adjustable pressure control to accommodate the patient's weight.