

Low Air Loss Alternating Pump & Mattress System
ALM831-842T



Emergency CPR Valve Located on Side

- Controls pressure, moisture, shear and friction, for prevention and treatment of stages I-IV Pressure sore ulcer
- Dual bladder design prevents patients from bottoming out while sitting up and will maintain air during power outages
- Adjustable anchor straps secures mattress to bed frame
- Multi-direction stretch fabric is zipper on, it is easy to remove for washing
- Quiet, dependable digital air controller unit provides hours of trouble free operation
- Alternating modes are adjustable in 4 timing selections, from 5, 10, 15, to 20 mins.
- Optional fast filling air blower
- HCPCS Code: E0277

Model	Overall Width	Overall Length	Overall Height	Air Cells	Power	Auto-Weight Detection	Weight Capacity
ALM831T	36"	80"	8"	13	110V / 60Hz	Yes	350 lbs.
ALM842T	42"	80"	8"	13	110V / 60Hz	Yes	600 lbs.

West: 800-367-6160

East: 800-548-6596